




proGAV[®] 2.0

IN TOUCH WITH YOU

proGAV[®] 2.0

VALVE AND SOFT TOUCH TOOLS



MIETHKE gravitational valves provide a reliable solution in prevention of overdrainage complications^{1,2} associated with the management of hydrocephalus.

proGAV 2.0 offers the reliability of our existing technology with a new tactile feedback mechanism and soft-touch tools.

proGAV 2.0 represents the latest development in refining our existing *proGAV* technology to better meet the needs of our customers and their patients.



The adjustable differential pressure unit is equipped with a new valve surface: By using slight pressure with your finger, the active lock mechanism of the valve is released and a tactile feedback is received.


The combination of the refined *proGAV 2.0* technology and soft-touch tools provides comfort for patients and security against inadvertent readjustment caused by magnet fields up to 3 Tesla.³

With the new *proGAV 2.0 Compass* the opening pressure of the valve can be identified quickly, without having to expose the patient to X-ray.

- 1 Thomale U.W., Gebert A.F., Haberl H., Schulz M. Shunt survival rates by using the adjustable differential pressure valve combined with a gravitational unit (*proGAV*) in pediatric neurosurgery. *Childs Nerv Syst.* 2013 Mar; 29(3):425-31.
- 2 Lemcke J., Meier U., Müller C., Fritsch M.J., Kehler U., Langer N., Kiefer M., Eymann R., Schuhmann M.U., Speil A., Weber F., Remenez V., Rohde V., Ludwig H.C., Stengel D. Safety and efficacy of gravitational shunt valves in patients with idiopathic normal pressure hydrocephalus: a pragmatic, randomised, open label, multicentre trial (SVASONA). *J Neurol Neurosurg Psychiatry.* 2013 Aug; 84(8):850-57.
- 3 Allin D.M., Czosnyka Z.H., Czosnyka M., Richards H.K., Pickard J.D. In vitro hydrodynamic properties of the Miethke *proGAV* hydrocephalus shunt. *Cerebrospinal Fluid Res.* 2006 Jun; 3:9. doi:10.1186/1743-8454-3-9.

proGAV[®] 2.0

FUNCTIONALITY OF VALVE AND BODY POSITION



The adjustable differential pressure unit can be set in steps of 1 cmH₂O between 0 and 20 cmH₂O. The gravitational unit is available in various pressure settings.

FUNCTIONALITY

The *proGAV 2.0* is a posture-dependent valve, i.e. its opening pressure changes gradually to correspond with the patients body position. For complete customization of the patients needs, the *proGAV 2.0* allows a specific opening pressure when the patient is lying down, and an opening pressure for when the patient is upright.


SUPINE POSITION

When the patient is in the supine position, only the adjustable differential pressure unit is active, which is preset to 5 cmH₂O*.

The gravitational unit is not active in this body position.



The functionality of the *proGAV 2.0* is illustrated interactively in the MIETHKE App.



UPRIGHT POSITION


In the upright position, the gravitational and adjustable differential pressure unit work together. As the patient gets in upright position, the tantalum ball (green) is activated within the gravitational unit and causes an increase in the valve opening pressure, due to its gravitational force.

The gradual increase in opening pressure in the upright position effectively protects against excess drainage.^{4,5}

In the example shown, a gravitational unit with 25 cmH₂O has been selected. The total opening pressure therefore amounts to 30 cmH₂O** while standing.

4 Lemcke J., Meier U., Müller C., Fritsch M., Kiefer M., Eymann R., Kehler U., Langer N., Schuhmann M.U., Speil A., Weber F., Remenez V., Rohde V., Ludwig H.C., Stengel D. On the method of a randomised comparison of programmable valves with and without gravitational units: the SVASONA study. *Acta Neurochir Suppl.* 2012; 114:243-46.

5 Freimann F.B., Sprung C. Shunting with gravitational valves - can adjustments end the era of revisions for overdrainage-related events?. *J. Neurosurg.* 2012 Dec; 117(6):1197-204.



After the adjustable differential pressure unit has been localized with the finger, the *proGAV2.0 Compass* is applied over the valve in the direction of flow.

The closed *proGAV2.0 Compass* indicates the opening pressure setting.

LOCALIZATION AND PRESSURE SETTING IDENTIFICATION

The *proGAV2.0 Compass* is used to localize the adjustable differential pressure unit and identify current setting of the valve.

The compass needle reflects the current opening pressure on the scale – Compass lid is closed.

The Compass should be aligned over the valve with the aid of the integrated template – For this purpose, the compass lid is opened.



The *proGAV 2.0 Adjustment Tool* centered over the valve.



By pressing lightly with the finger on the adjustable differential pressure unit, the mechanical “active lock” mechanism is deactivated and the opening pressure is set.

Releasing the finger pressure automatically locks the valve.


ADJUSTING THE VALVE

With the *proGAV 2.0 Adjustment Tool*, the opening pressure of the adjustable differential pressure unit can be set from 0 to 20 cmH₂O.

The *proGAV 2.0 Adjustment Tool* has to be placed above the valve, so that the value of the desired opening pressure points towards the valve inlet connector.

proGAV[®] 2.0

+ proGAV 2.0 valve with SA 2.0



Connector: $d_o = 1.9$ mm

Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX642T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX643T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old


Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
FX410T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX640T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX641T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX644T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX645T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0 with Distal Catheter

+ proGAV 2.0 valve with Distal Catheter, 1200 mm



Connector: $d_o = 1.9$ mm

Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX648T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX649T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
FX417T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX646T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX647T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX650T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX651T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O


** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0



Shunt System with Pediatric CONTROL RESERVOIR

- + proGAV 2.0 valve with Pediatric CONTROL RESERVOIR* and Distal Catheter, 1200 mm
- + Ventricular Catheter with Pediatric Burrhole Deflector and introducing stylet, 250 mm

* An additional valve in the base of the Pediatric CONTROL RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the Ventricular Catheter and the distal drainage.



Pediatric
CONTROL RESERVOIR



Connector: $d_o = 1.9$ mm

Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX609T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX610T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available


Art. no.	Adjustable unit	Gravitational unit
FX606T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX607T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX608T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX611T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX612T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

Shunt System with CONTROL RESERVOIR

- + proGAV 2.0 valve with CONTROL RESERVOIR* and Distal Catheter, 1200 mm
- + Ventricular Catheter with Burrhole Deflector and introducing stylet, 250 mm

* An additional valve in the base of the CONTROL RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the Ventricular Catheter and the distal drainage.



Connector: $d_o = 1.9$ mm
 Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX602T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX603T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
FX431T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX600T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX601T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX604T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX605T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0



Shunt System with Pediatric SPRUNG RESERVOIR

- + proGAV 2.0 valve with Distal Catheter, 1200 mm
- + Pediatric SPRUNG RESERVOIR* with Distal Catheter, 600 mm
- + Ventricular Catheter with introducing stylet, 180 mm

* An additional valve in the base of the Pediatric SPRUNG RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the Ventricular Catheter and the distal drainage.



Pediatric
SPRUNG RESERVOIR



Connector: $d_o = 1.9 \text{ mm}$

Catheter: $d_i = 1.2 \text{ mm}$, $d_o = 2.5 \text{ mm}$

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX583T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX584T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
FX580T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX581T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX582T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX585T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX586T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O


** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0

Shunt System with SPRUNG RESERVOIR

- + proGAV 2.0 valve with Distal Catheter, 1200 mm
- + SPRUNG RESERVOIR* with Distal Catheter, 600 mm
- + Ventricular Catheter with introducing stylet, 180 mm

* An additional valve in the base of the SPRUNG RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the Ventricular Catheter and the distal drainage.



Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX576T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX577T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
FX424T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX574T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX575T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX578T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX579T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0



Shunt System with Pediatric SPRUNG RESERVOIR

- + proGAV 2.0 valve with integrated Pediatric SPRUNG RESERVOIR* and Distal Catheter, 1200 mm
- + Ventricular Catheter with introducing stylet, 180 mm

* An additional valve in the base of the Pediatric SPRUNG RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the ventricular catheter and the distal drainage.



Pediatric
SPRUNG RESERVOIR



Connector: $d_o = 1.9$ mm

Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX636T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX637T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available


Art. no.	Adjustable unit	Gravitational unit
FX633T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX634T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX635T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX638T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX639T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

Shunt System with SPRUNG RESERVOIR

- + proGAV 2.0 valve with integrated SPRUNG RESERVOIR* and Distal Catheter, 1200 mm
- + Ventricular Catheter with introducing stylet, 180 mm

* An additional valve in the base of the SPRUNG RESERVOIR makes it possible to flush the fluid only in the distal direction. This feature allows for controlling the patency of the ventricular catheter and the distal drainage.



Connector: $d_o = 1.9 \text{ mm}$
 Catheter: $d_i = 1.2 \text{ mm}$, $d_o = 2.5 \text{ mm}$

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX629T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX630T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available


Art. no.	Adjustable unit	Gravitational unit
FX626T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX627T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX628T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX631T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX632T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0

Shunt System with Pediatric Burrhole Reservoir

- + proGAV 2.0 valve with Distal Catheter, 1200 mm
- + Pediatric Burrhole Reservoir with Distal Catheter, 600 mm
- + Ventricular Catheter with introducing stylet, 180 mm



Connector: $d_o = 1.9$ mm

Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
FX570T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
FX571T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available



Art. no.	Adjustable unit	Gravitational unit
FX445T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
FX568T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
FX569T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
FX572T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
FX573T	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0

Shunt System

- + proGAV 2.0 valve with Distal Catheter, 1200 mm
- + Ventricular Catheter with Pediatric Burrhole Deflector and introducing stylet



Connector: $d_o = 1.9$ mm
 Catheter: $d_i = 1.2$ mm, $d_o = 2.5$ mm

Standard configurations

Art. no.	Adjustable unit	Gravitational unit	Recommendation**
8220 2020	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	20 cmH ₂ O	Children up to 5 and adults over 65 years old
8220 2025	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	25 cmH ₂ O	Individuals between 5 and 65 years old

Alternative configurations available

Art. no.	Adjustable unit	Gravitational unit
8220 2000	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	without
8220 2010	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	10 cmH ₂ O
8220 2015	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	15 cmH ₂ O
8220 2030	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	30 cmH ₂ O
8220 2035	0 - 20 cmH ₂ O (preset to 5 cmH ₂ O)	35 cmH ₂ O

** Recommendations should only be used as guideline as other settings may also be necessary depending on the patient and anatomy.

proGAV[®] 2.0 Tools

+ proGAV 2.0 Adjustment Tool, proGAV 2.0 Compass, proGAV 2.0 Tool Set and proGAV 2.0 Checkmate



proGAV 2.0 Adjustment Tool



proGAV 2.0 Compass



proGAV 2.0 Tool Set



proGAV 2.0 Checkmate

Art. no.	Tools
FX400T	<i>proGAV 2.0 Adjustment Tool</i>
FX401T	<i>proGAV 2.0 Compass</i>
FX404T	<i>proGAV 2.0 Tool Set (contains FX400T and FX401T)</i>
FX409T	<i>proGAV 2.0 Checkmate</i>
on request	<i>proGAV 2.0 X-ray template and pressure recommendation card</i>



Manufacturer:

Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Tel. +49 331 62 083-0 | Fax +49 331 62 083-40 | www.miethke.com

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.
Detailed information can be requested from info@miethke.com.
Technical changes reserved. This brochure may be used exclusively for providing information about our products. Reprinting, even in excerpts, is prohibited.