

# Capre MC Medical Chair

Suitable for dialysis and cytostatic treatments

Lojer Capre MC Medical Chair improves patient comfort during various long-term treatment. Versatile adjustments and number of quality accessories facilitate the positioning of the patient.



The luxurious Capre MC multipurpose chair has been designed together with healthcare professionals to fulfil the demands of today's healthcare processes in procedures like dialysis and cytostatic treatment. The Capre MC chair features pleasantly contoured, soft imitation leather, individual adjustment possibilities and built-in comforts, such as an iPad stand, tray, IV-pole and reading light. All of this combines to improve patients' quality of life during treatment, whilst also improving the quality image of the healthcare facility. The control system allows the independent adjustment of the chair height, back section, seat angle and leg sections. The leg section is also motorized. The head section also feature a length adjustment mechanism. The chair features five pre-programmed positions: three treatment positions (sitting, beach chair and lying), Trendelenburg position and one memory slot for storing the customer's own settings. Due to its modular design, the chair can also be used as an examination chair.



Scan the QR code  
to see the product video, or visit  
<http://youtu.be/EbrNPcGUhyk>

**LOJER**® *For easy care*

## Technical specifications

- Total length: 186..213 cm
- Height adjustment: 60..85 cm
- Castors: Ø 75mm, centrally locking
- Back section adjustment: -12°..+90°
- Head section adjustment: 0°..10°
- Back section length: 600 mm
- Seat adjustment: +12°..0°
- Seat width: 620 mm
- Total width: 880 mm
- Leg section adjustment: -90°...0°
- Arm support adjustment: 0°..180° rotation
- Safe Working Load (SWL): 160 kg
- USB output: Power 4W
- Electric system: 24 V
- CE marking Yes
- Warranty: 2 years
- Country of Manufacture: Finland



Features	
CE -marked medical device	x
Luxurious anatomically contoured upholstery	x
Two colour options for the plastic parts	x
Upholstery SKAI Fashion (16 options)	x
Design castors for moving the chair, centrally lockable	x
Intelligent control system with pre-programmed positions	x
Trendelenburg (12° max.)	x
Motorised adjustment of back section	x
Motorised adjustment of seat angle	x
Motorised height adjustment	x
Motorised one-piece leg section with optional foot support	x
Head section length and angle adjustment	x
Arm rests stay horizontal regardless of back section angle	x
180° rotation of arm rests	x
Safe, quiet 24V electrical system, USB charger	x
Warranty 2 years	x
Optional features:	
Versatile arm support with angle and height adjustment – available in two different cushion size (standard and large with concave shape)	o
Foot support with stepless hydraulic length adjustment	o
Coffee tray	o
Tablet PC stand with flexible neck	o
Face cut-out	o
Plug for face cut-out	o
Decorative stitching in the upholstery	o
LED light (usb)	o
Infusion stand	o
Wired hand control (instead of flexible arm)	o
Foot controller	o

Explanation of markings: x = standard feature, o = optional accessory, - = not available

## Accessories



Coffee tray



Tablet support with flexible arm



Versatile arm support with large, concave cushion



Led light



Foot controller

**LOJER**® For easy care

Lojer Oy, Putajantie 42, P.O. Box 54, FI-38210 Sastamala, Finland  
Tel. +358 (0)10 830 6700 | info@lojer.com | www.lojer.com