

Delta Treatment Tables

Excellent Value for Money

This classic massage table has been the favourite of masseurs for a long time. Delta Professional massage tables represent the best in their field, both in terms of function and durability. Our many years of experience can be seen in the design of our tables, ensuring stability and durability at all working heights and positions.

Professional and Standard. Standard treatment tables are basic tables for professionals; more affordable than the Professional series, nevertheless with features that fulfil the needs of demanding users.

"Functional, stable and durable."

Delta Standard DS5H Treatment table with castors



Delta Professional DP4 Massage table with castors, pair of side supports and ankle cushion

LOJER[®] For easy care

Delta Professional	DP2	DP3	DP 4	DP5				
Sections	2	3	4	5				
Width, 50, 55 or 60 cm	X	X	X	X				
Width, 70 cm (side supports cannot be used)	X	-	X	-				
Length, 195 cm	X	X	X	X				
Weight, kg	75	76	76	78				
Height adjustments, cm	47 - 98	48 - 100	47 - 98	48 - 100				
Hydraulic height adjustment	-	-	-	-				
Safe working load (SWL), 200kg	X	X	X	X				
Foot section adjustment range	-	0 - 75°	-	0 - 75°				
Middle section adjustment range	-	60°	-	60°				
Electric drainage and foot section adjustment	-	•	-	•				
Centrally lockable castors	•	•	•	•				
Side supports	•	•	•	•				
360° height adjustment foot bar	•	•	•	•				
Double-hinged face section	•	•	•	•				
Hand controller	•	•	•	•				
Foot controller	•	•	•	•				
Foot switch support	•	•	•	•				
Paper roll holder	•	•	•	•				
Safety switch to prevent accidental use	•	•	•	•				
SKAI Fashion upholstery (13 options)	X	X	X	X				
Electrical system, 230 V	X	X	X	X				
CE-marking	X	X	X	X				
Warranty 2 years	X	X	X	X				
Country of Origin	Finland	Finland	Finland	Finland				
Delta Standard	DS2	DS2H	DS3	DS3H	DS4	DS4H	DS5	DS5H
Sections	2	2	3	3	4	4	5	5
Width, 50, 55 or 60 cm	X	X	X	X	X	X	X	X
Width, 70 cm (side supports cannot be used)	X	X	-	-	X	X	-	-
Length, 195 cm	X	X	X	X	X	X	X	X
Weight, kg	60	60	61	61	62	62	64	64
Height adjustments, cm	57 - 87	57 - 87	58 - 89	58 - 89	57 - 87	57 - 87	58 - 89	58 - 89
Hydraulic height adjustment	-	X	-	X	-	X	-	X
Safe working load (SWL), 200kg	X	X	X	X	X	X	X	X
Foot section adjustment range	-	-	0 - 70°	0 - 70°	-	-	0 - 70°	0 - 70°
Middle section adjustment range	-	-	60°	60°	-	-	60°	60°
Electric drainage and foot section adjustment	-	-	•	•	-	-	•	•
Centrally lockable castors	•	•	•	•	•	•	•	•
Side supports	•	•	•	•	•	•	•	•
360° height adjustment foot bar	•	-	•	-	•	-	•	-
Double-hinged face section	•	•	•	•	•	•	•	•
Hand controller	•	-	•	-	•	-	•	-
Foot controller	•	-	•	-	•	-	•	-
Foot switch support	•	-	•	-	•	-	•	-
Paper roll holder	•	•	•	•	•	•	•	•
Safety switch to prevent accidental use	•	-	•	-	•	-	•	-
SKAI Fashion upholstery (13 options)	X	X	X	X	X	X	X	X
Electrical system, 230 V	X	-	X	-	X	-	X	-
CE-marking	X	X	X	X	X	X	X	X
Warranty 2 years	X	X	X	X	X	X	X	X
Country of Origin	Finland	Finland	Finland	Finland	Finland	Finland	Finland	Finland

X = Standard - = not available • = Optional



360° height adjustment foot bar



Double-hinged face section



Side supports