

 CRISTALENS

LUCIS[®]

HYDROPHILIC LENS



► Monofocal

MADE
IN FRANCE

CE 0459

TECHNICAL SPECIFICATIONS

Lens type	For implantation in the capsular bag
Optic diameter	6.00 mm
Overall diameter	11.00 mm (from -9.0D to +9.0D) 10.50 mm (from +10.0D to +30.0D)
Design	One-piece
Optic design	Monofocal Spherical
Angulation	10°
Material	Hydrophilic Benz 25 made of hydrophilic and hydrophobic copolymers
Dioptic powers	From -9.0D to +10.0D by 1.0D From +10.0D to +30.0D by 0.5D
Estimated A-Constant (SRK-T)	118.0 Ultrasound biometry 118.5 Interference laser biometry
Suggested anterior chamber depth (ACD)	4.96 mm Ultrasound biometry 5.25 mm Interference laser biometry
Refractive index	1.46
Sterilization	Steam
Recommended injection system	Single-use injector
Recommended incision size	2.2 mm



www.cristalens-international.com

 **CRISTALENS INDUSTRIE**

4 rue Louis de Broglie - 22300 Lannion - FRANCE
Tel +33 (0)2 96 48 92 92 - Fax +33 (0)2 96 48 97 87

