

Reanibex 100

Feel Cardioprotected

2 options:

- Semiautomatic
- Fully- Automatic

Adult and pediatric patient use

Alerts about the battery and pads status

CPR with metronome and feedback

Sigfox Transmission

Geolocation

Configured according to ERC/AHA guidelines

Reanibex 100

Ready for life



FEATURES



CLEAR GUIDANCE

Audible prompts to guide users through the complete operation sequence. Volume adapted to background noise.



WITH COMMUNICATIONS

Equipment and accessories status alerts transmission via Sigfox and Geolocation.



SEMI- OR FULLY-AUTOMATIC

In the Semi-automatic version the rescuer has to push the shock button to deliver the shock. In the fully automatic version the shock is delivered automatically when needed.



FASTER SHOCK

Automatically switched on when the cover is opened, pre-connected pads and precharge for a faster shock delivery.



SELF CHECK

Automatic daily selftest to ensure device and accesories continuous readiness.



CPR SUPPORT

Metronome and rate feedback in real time for effective chest compressions.

APPLICATIONS



REANIBEX DATA NET

Web application for monitoring AED devices. Alerts generation about the status of the devices and their accessories, equipment configuration, software updates and the possibility to display test results and episodes (episodes optional).



REANIBEX NFC READER

NFC application for Android devices to check the status of the REANIBEX 100 and its accessories without PC connection.



REANIBEX USB SYNC

PC application to synchronize the AED with the Reanibex Data Net via USB to configure the device and to download the information stored in the REANIBEX 100.

Bexen cardio

Osatu S. Coop. Edificio Zearrekobuelta Subida Areitio, 5 48260 Ermua (Spain) Tel: +34 943 170 220 Fax: + 34 943 170 227 info@bexencardio.com www.bexencardio.com

