

Bexen cardio

Ready for life

Reanibex 300 MANUAL

AED with Manual override

Biphasic with Manual override
up to 360J

Universal AED algorithm for
adult and paediatric patients

CPR guidance device in real time
(Push Pad)

Possibility to configure the
screen on high contrast mode

“Autotraining” mode in the same
equipment



Reanibex 300 MANUAL.

User friendly and very intuitive.



FEATURES

- Biphasic technology.
- Universal AED algorithm for adult and paediatric patients.
- Large colour display 5.7".
- Possibility of setting the screen on high contrast mode.
- "Push Pad" device, reports the quality of the frequency and depth of chest compressions in real time (optional).
- Autotraining mode (optional).
- Configurable CPR protocol.
- Configured according to ERC/ AHA guidelines 2015.
- Automatic daily selftest ensures the device continual readiness.
- USB connection.
- Rugged and durable design for any environment.

ACCESSORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Adult and paediatric ECG electrodes.
- Carrying case (1).
- Ambulance bracket (EN 1789:2007 certified) (2).
- 3 leads patient cable.
- "Push Pad" device.
- Rechargeable and non rechargeable battery (3).
- External battery charger.
- External dual battery charger 12 V (DC) or 220V (AC).
- "REANIBEX Data Manager" Software application to manage and organize all the data gathered.
- Rechargeable battery for autotraining mode (4).



Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobuelta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com