

Bexen cardio
Ready for life

Reanibex 800

The unique modular
Defibrillator Monitor
“plug & play”
of the market



Defibrillator Monitor
multiparameter

Manual and AED up to 360J

Universal AED algorithm for
adult and paediatric patients

12 leads acquisition

ECG interpretation

Non invasive pacing

SpO2, SpCO, SpMet, SpHb,
SpOC, PVI

NIBP

EtCO2

IBP

TEMP

CPR assistance

Data transmission
(Bluetooth, GSM)

Reanibex 800. Defibrillator Monitor modular.



Communications (Reanibex Data Manager, Data Link, Data Cloud)



Operates from the mains (AC), from a vehicle battery (DC) and with its battery

12 ECG acquisition key



Modular easily upgradeable. Plug & Play system

50 mm printer

Non invasive pacing

External reusable paddles with contact indicator on the APEX paddle



Status indicator

Large colour display (8.4")

Therapy selector

AED for adult or paediatric

FEATURES

- Biphasic technology.
- Large colour display (8.4").
- Monitoring and diagnostic capabilities.
- 3, 5 or 10 leads patient cable.
- Non invasive pacing (optional).
- 12 leads ECG data transmission from the field to the Emergency Service (optional).
- 12 leads Interpretive Algorithm (University of Glasgow) (Optional).
- Synchronized Cardioversion.
- CPR Metronome with feedback.
- Capnography from Oridion with Microstream technology (Optional).
- Non-invasive Blood Pressure from SunTech Medical.
- Masimo Rainbow SET technology that monitors SpO₂, SpCO and SpMet, SpHb, SpOC, PVI (optional).
- Invasive Blood Pressure (two lines).
- Temperature (two lines).
- Operates from the mains (AC), from a vehicle battery (DC) and with its battery.
- Possibility to operate with two batteries.
- Possibility of setting the screen on high contrast mode.
- Transmission (GSM, Bluetooth) to Emergency Services, Hospitals... configurable according to the user needs (optional).

ACCESSORIES

- Adult and paediatric disposable pads for defibrillation.
- Adult and paediatric ECG electrodes.
- Carrying case (1).
- Neoprene upper bag.
- Disposable internal paddles (The iD-System) (2).
- 3, 5 or 10 leads patient cable.
- Rechargeable battery.
- External battery charger.
- Connection cable at 12V (DC).
- Ambulance bracket (EN 1789:2007) (3).

Applications

- REANIBEX Data Manager, Data Link and Data Cloud software applications to manage and organize all the data (4).




(1)



(3)



(2)



(4)

Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobuelta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com

