



1) Rapid Warm-up & Dry System




2) Convenient and Safe Sterilizing Agent



3) 7" full touch LCD, Easy Monitoring the Cycle Information



4) Plug and Play



5) All-in-one Basket Tray



Sterilizing Agent		Hydrogen Peroxide 25 Cycles/Bottle
Total Cycle Time	Quick	25+5min
	Standard	35+5min
	Special	45+5min
Cycle Temperature		50 ± 5°C
SAL(Sterility Assurance Level)		10 ⁻⁶
By-products		Oxygen and vapor water only
Lumenclaims		Rigid lumen & Flexible lumen, Endoscopes
Chamber	Type Material	Circle Stainless Steel (SUS)
Dimensions	Overall Chamber	680mm(W) X 1575(H) X 900(D) Ø420mm X 720(D)
Volume		Total : 100 Liters
Weight		250Kg
Control		Microprocessor & Windows CE embedded
Cycle Information		Screen, Printer, USB, 100/10Mbps Ethernet(Optional)
PM Cycle		Automatic Alarm & Setup
Electrical		110V / 220V, 50 / 60Hz, 1Phase, 3000W
Installation Requirements		Front : 100cm Rear, Left side, Right side : 10cm Placement : Built-in wheels provide mobility
Room Conditions		5~40°C, 0~95%RH(Non-Condensing)
Printer		Built-in Thermal Printer(60Ø Roll Paper) Cycle Parameters(Temp, Pressure, Time, Daily&Total Cycle, etc) Alarm & Error Display
Others		Emergency Stop(Front), Operator ID Login, Self-Diagnosis, Multi-Language, RFID Coding System