





## 1) Rapid Warm-up & Dry System




## 2) Convenient and Safe Sterilizing Agent



## 3) 7" full touch LCD, Easy Monitoring the Cycle Information



## 4) Plug and Play



## 5) Table-top Model



## 6) Oilless Pump



## 7) Useful Cart



Sterilizing Agent		Hydrogen Peroxide 40 Cycles/Bottle
Total Cycle Time	Quick	25+5min
	Standard	35+5min
	Special	45+5min
Cycle Temperature		50 ± 5°C
SAL (Sterility Assurance Level)		10 <sup>-6</sup>
By-products		Oxygen and vapor water only
Lumenclaims		Rigid lumen & Flexible lumen, Endoscopes
Chamber	Type Material	Rectangular Stainless Steel (SUS)
Dimensions	Overall Chamber	<b>530mm(W) X 500(H) X 735(D)</b> <b>260mm(W) X 220(H) X 590(D)</b>
Volume		<b>Total : 33 Liters</b>
Weight		<b>90 Kg</b>
Control		Microprocessor & Windows CE embedded
Cycle Information		Screen, USB, 100/10Mbps Ethernet (Option)
PM Cycle		Automatic Alarm & Setup
Electrical		<b>110V / 220V, 50 / 60Hz, 1Phase, 1500W</b>
Installation Requirements		Front : 100cm Rear, Left side, Right side : 10cm Placement : Built-in wheels provide mobility
Room Conditions		5-40°C, 0-95%RH (Non-Condensing)
Others		<b>Operator ID Login, Self-Diagnosis, Multi-Language, RFID Coding System</b>
Cart (Option)		<b>530mm(W) X 760(H) X 700(L) / 30Kg</b>