



NEW


# SONIC SHOT<sup>®</sup> 7

MR CONTRAST DELIVERY SYSTEMS

MR Medium Contrast Injector Sonic Shot 7

# Sonic Shot 7 Sharpened Functionality

Each performance specification was honed to perfection, for safe and reliable contrast study.



## ■ Ultrasonic motor equipped

Completely non-magnetic design. Exclusive for MR, long-lifetime and trusting performance

## ■ Powerhead Lamps

Quickly & Easily recognize injector operation



## ■ Auto Detect Plunger

Insert the syringe, and the injector will automatically detect the plunger position. No need to manipulate manually.



## ■ One-Touch Syringe Adaptor

Exclusive adaptors designed for pre-filled syringes. Easily inserted and removed with One-Touch design



## ■ Operate Button

Powerhead buttons are designed for easy identification and ease of operation.

## Power Head

Completely non-magnetic structure, including new features, and enhanced ease of use

■ Body Weight Input Mode



Injector will automatically calculate the protocol based only on the patients weight

■ Body Weight Mode




■ Drip Infusion



Simplified Drip Mode, easy access and easy to setup and use

■ injection



■ Scan Timing



Real-time information, Elapsed Time bargraph and audible sound indicate when to scan



# Console

Design for improved reliability and safety with Nemoto's advanced Body Weight Protocol.

