


High definition, high efficiency

Optimum combination of high performance and multi-function.



Large size flat panel detector brings clear image with rich details

Precision imaging with no distortion

Self-developed high voltage generator

1.0ms


Minimum exposure time

1000mA

Maximum current output

80kW

High voltage generator power



Diagnostic table, smart control

Multi functionality and outstanding technics

220cm

Tabletop length

120cm

Longitudinal movement

100cm

Minimum SID

Easy operation

Complex moving controlled in one hand

Unique design of control handle, brings an integration of table moving control and spot film acquisition control. Buttons function can be defined by the operator according to his preference.



Intelligent control handle



One handle to control all movements



Intelligent image stitching

Intelligent and auto stitching of full spine, arm or legs





Dynamic / Static image stitching



Auto program controlled stitching



Fast and precision imaging



Multifunctional machine

Digital fluoroscopy, spot image, DR

Artistic design, easy operation

