



Stepping into 100μm vision

- 100μm large size detector brings precise imaging
- Ten-axis integrated intelligent automatic positioning
- Six-way intelligent tracking positioning system
- InvaRay intelligent imaging platform

80kW 350kHU 17x17"

100μm large size detector / Precise imaging



- 2.4/5G Dual Frequency Wireless Transmission
- 4.3K×4.3K (18.5M pixels)
- 8s imaging
- 43×43cm
- DQE 73%

4K×4K, 18.5M pixel detector / improves the detail detection capacity



Ten-axis integrated intelligent automatic positioning / Six-way intelligent tracking system

- # Intelligent cross-rail suspended tube stand
- ⊞ Ten-axis integrated intelligent automatic positioning
- Ⓜ Elevating RAD table
- ⊕ Four-way floating tabletop

Auto positioning / Auto tracking / Auto sync / Auto stitching
Quick switch to pre-programmed positions




Intelligent touch system / Easy operation, powerful function



- Patient information display
- Exposure parameter setting
- Last image preview
- Safety warning
- APR selection
- Patient positioning illustration
- Historical image retrieve/review
- Gravity adaptive

Intelligent one touch automatic stitching / Accurate and convenient



Dual WiFi Detectors, flexible in different practices / Whole body radiography without detector exchange, enhance work flow



Intelligent diagnosis / Integrate artificial intelligence to quickly identify lesions and assist doctors in accurate diagnosis



Intelligent bone suppression function

Intelligent knee joint diagnosis

Intelligent diagnosis of tuberculosis

Intelligent pneumothorax diagnosis

Intelligent image search function