





Mobile radiography-Ingenious design and long battery life

-  Suit for radiological department, ward, ICU and operating room
-  High power, competent for different clinical demands
-  Large capacity battery pack, to ensure long battery life
-  Wireless large size detector, high definition image acquisition at low dose




Strong Power

Comparable to the stationary X-ray machine

40kW 500mA 125kV 400mAs

Outstanding imaging performance




Presents good quality image at low dose

-  Wireless large size detector, easy to use
-  High gray scale, better image contrast
-  Optional sizes: 35cmx43cm, 43cmx43cm



Long battery life

Fully charged battery can support more than one thousand times of exposures

-  Large capacity battery, long standby
-  Quick charge, suit for heavy duty work
-  Intelligent battery management, extend the life of battery

Smart and compact

Slim in shape, more flexible moving




250 kg 1470mm 610 mm
Weight height width



Artistic design,
easy operation

One touch operation

Simple and convenient

-  Hold to move, thanks to the tactile switch of wheels.
-  Detachable Surface pro, provides convenient touch-screen control.
-  Wireless hand switch, for remote exposure control.

