



Intimate care


Digital mammography solution

- Flat panel detector, full coverage, high efficiency examination
- Intelligent electric compression device, comfort and safety
- Advanced imaging algorithm, low radiation dose, high definition image

Dedicated flat panel detector for mammography

Small lesion, clearly shown

- Fast imaging
- Large imaging area
- Low radiation
- Magnification radiography



Dedicated X-ray tube for mammography

Low radiation dose, high definition image

- Suitable for breast examine
- Imported Molybdenum-target X-ray tube
- Long life X-ray tube
- Molybdenum and Rhodium filtration, auto-switching

Integrated design of mechanical structure

--- stable, flexible and easy to operate

- Angle memory function, fast and precise positioning
- Concentric rotation
- One key switching to programmed position
- 270 degree rotation, meet the clinic needs



Intelligent electric compression device

Gentle care during the whole procedure of examination

- Compression strength will automatically adjust according to breast tissue density
- Support manual compressor release
- Auto memory of compression strength and angle, enable positioning fast and precise

Smart operation program

- One switch moving to the programmed position
- AEC exposure control
- Foot switch controlled compressor movements, more convenient and more efficient



Intelligent software platform

Intelligent control brings ease of use and high efficiency

- Advanced imaging algorithm
- Modular configuration
- DICOM 3.0
- Friendly acquisition interface
- Smart configuration of acquisition protocol