



Optimum integration of top-level hardware

Self-developed core components

Be able to supply sustainable service and system upgrade



Top quality magnet

Patent design


Using highest level magnetic materials

- Reliable open structure
- Unique nano anti-eddy current technology
- Stable magnetic field, no need of manual shimming


Digital optical fiber spectrometer

Improves SNR significantly

Magnet room



Magnet room



RF receiving coils

Specialized coils can be customized according to customer demands

- Phase array coil** High acquisition speed High definition image
- Flexible coil** Care for obese patient
- Customized coil** Customized coils can be developed for specialized demands



Intelligent imaging platform


Dual screen workstation

- Windows based interface
- Multi-language support
- 9 scout images pilot scan
- Classified management for user-defined sequence library
- Multi-task quick switching
- Corner shortcuts for quick image processing
- Image navigation view

Pad review and control

Anytime, anywhere

- Image auto-sending to the PACS
- Convenient for remote image review



Cloud diagnosis service (for contracted customers only)

- Wanli Cloud diagnosis service provides professional image report and extended services
- Remote image review and data log reading
- Remote service system, provides real-time equipment status monitoring