


## Optimum integration of top-level hardware

### Self-developed core components


Be able to supply sustainable service and system upgrade



## Open magnet system

### Patent design


- Top-level magnetic material
- Small size, light weight
- 320° wide open C-shape design
- Constant temperature technology



## Digital optical fiber spectrometer

### Improves SNR significantly

Magnet room



## RF receiving coils

### Specialized coils can be customized according to customer demands

- Phase array coil** High acquisition speed High definition image
- Flexible coil** Care for obese patient
- Customized coil** Customized coils can be developed for specialized demands



## Intelligent imaging platform

### Dual screen workstation

- Windows based interface
- Multi-language support
- 9 scout images pilot scan
- Classified management for user-defined sequence library
- Multi-task quick switching
- Corner shortcuts for quick image processing
- Image navigation view



## Pad review and control

### Anytime, anywhere

- Image auto-sending to the PACS
- Convenient for remote image review



## Cloud diagnosis service (for contracted customers only)

- Wanli Cloud diagnosis service provides professional image report and extended services
- Remote image review and data log reading
- Remote service system, provides real-time equipment status monitoring

